

ANSWER


1. Which is the correct reading of the meniscus as shown in the portion of the buret below?



- (1) 41.30 ml (2) 41.35 ml
- (3) 42.60 ml
- (4) 42.65 ml
- 2. Which device should be used to accurately measure a volume of 16.30 milliliters?










- 3. Which is a characteristic of an aqueous solution of HNO3? ACID
  - (1) It conducts electricity.
  - (2) It forms OH ions.
  - (3) It turns litmus blue.
  - (4) It turns phenolphthalein pink.
- 4. Which of the following is the best conductor of electricity?
  - (1) NaCl(s)
- $(3) C_6H_{12}O_6(s)$
- (4)  $C_6H_{12}O_6(aq)$

COVALENT

5. Beakers A, B, C, and D shown below each contain a different solution.



The bulb will glow when the conductivity apparatus is placed into which beakers?

- (1) A and B
- 3) A and D
- (2) B and C
- (4) Cand D
- 6. A student was given four unknown solutions. Each solution was checked for conductivity and tested with phenolphthalein. The results are shown in the data table below.

Solution	Conductivity	Color with Phenolphthalein	
Α	Good	Colorless	
В	Poor	Colorless	
(0)	Good	Pink BA	
D	Poor	Pink	

Based on the data table, which unknown solution could be 0.1 M NaOH? NaOH

(1) A(2) B

- 7. According to the Arrhenius theory, the acidic property of an aqueous solution is due to an excess of
  - (1) H<sub>2</sub>

(3) H<sub>2</sub>O

- (4) OH-
- 8. According to the Arrhenius theory, when a base dissolves in water it produces
  - (1) H<sup>+</sup> as the only positive ion in solution
  - (2) NH<sub>4</sub><sup>+</sup> as the only positive ion in solution
  - 3 OH as the only negative ion in solution (4) CO<sub>3</sub><sup>2-</sup> as the only negative ion in solution

NaOH > OH

S	What are the relative ion concentrations in an acid solution?  19. Which statement describes the characteristics of					cteristics of an	
	more H <sup>+</sup> ions than OH <sup>-</sup> ions (2) fewer H <sup>+</sup> ions than OH <sup>-</sup> ions		Arrhenius base?				
			(1) It changes blue litmus to red and has a pH less than 7.				
	<ul> <li>(3) an equal number of H<sup>+</sup> ions and OH<sup>-</sup> ions</li> <li>(4) H<sup>+</sup> ions but no OH<sup>-</sup> ions</li> </ul>			(2) It changes blue litmus to red and has a pH greater than 7.			
According to the Arrhenius theory of acids, citric acid in oranges and acetic acid in vinegar are classified as				(3) It changes red		has a pH less	
			than 7				
	acids because their aqui	eous solutions contain	(	greater than 7.			
	1) hydrogen ions H+				171 BCE 11		
	(2) hydrogen atoms	(4) hydroxide atoms	20.	The table below wa aboratory data.	s compiled from ex	xperimental	
11	A sample of Ca(OH) <sub>2</sub> is	considered to be an		aboratory data.			
	Arrhenius base because it dissolves in water to yield (1) Ca <sup>2+</sup> ions as the only positive ions in solution					pH RANGE	
	(1) Ca <sup><math>+</math></sup> ions as the only	y positive ions in solution y positive ions in solution				AT WHICH	
	(3) OH ions as the only	negative ions in solution		INDICATOR	CHANGE	CHANGE OCCURS	
	(4) H <sup>-</sup> ions as the only r	egative ions in solution		Bromthymol Blue	(yellow) - blue	6.2 - 7.6	
12	Which pH value indicate	s the most basic solution?		Thymol Blue	red — (yellow).	1.2 - 2.8	
	(1) 7	(3) 3		Methyl Orange	red - vellow	3.1 - 4.4	
	(2) 8	4 11			The Committee	7	
13.	An acidic solution could	nave a pH of		At what pH would a	Il throe indicators o	VELLO	u
	(1) 7	(3) 3	у	rellow?	in three indicators a	ppear as [	,
	(2) 10	(4) 14	(	1) 1.9	(3) 4.7	2.8 4.4	0.
14.	Adding 0.1 M NaOH to a	0.1 M solution of HCI will	(	2) 2.9	(4) 8.7		
	cause the pH of the solution	tion to	21. (	One sample of a so	lution with a nH of	10 is tested	
	(1) decrease	(3) remain the same	V	vith phenolphthaleir	and another sam	ple of this	
	(2) increase		S	olution is tested with	th litmus. In this so	lution the color	
15.	15. Which relationship is present in a solution that has a			1) blue and the ph	enolphthalein is pir	nk TABLE	M
	pH of 7?  (1) $[H^+] = [OH^-] = \{1, 0\}$ (3) $[H^+] < [OH^-]$		(	<ol><li>red and the phe</li></ol>	nolphthalein is pinl	<	,,,
			<ul><li>(3) blue and the phenolphthalein is colorless</li><li>(4) red and the phenolphthalein is colorless</li></ul>				
	(2) [H <sup>+</sup> ] > [OH <sup>-</sup> ]	(4) $[H^{+}] + [OH^{-}] = 7$	(4	4) red and the phe	nolphthalein is cold	orless	
16.	A solution at 25°C with a		22. A	n aqueous solution	turns litmus red.	The pH of this	
	(1) more H <sub>3</sub> O <sup>+</sup> ions than	OH <sup>-</sup> ions		olution could be			
	(2) fewer H <sub>3</sub> O <sup>+</sup> ions than an equal number of H	OHTIONS		1) 14 2) 11	(3) 8	TABLE M	
	(4) no H <sub>3</sub> O <sup>+</sup> ions or OH <sup>-</sup>	ions			(4))5		
17	As an aqueous solution b	2 A20		Which 0.1 M solution			
OH	As an aqueous solution by hydroxide ion concentration	on		1) HBr(aq) Aco	(3) LiOH	(aq) BASE	
•	decreases	(3) remains the same	(2	2) CO <sub>2</sub> (aq)	(4) CH <sub>3</sub> C	H(aq) ALCOHUL	
	(2) increases		24. A	n indicator was use	ed to test a water s	olution with a	
18.	When the pH of a solution	changes from a nH of 5 to		H of 12. Which indi  I) colorless with lite		be observed?	
18. When the pH of a solution changes from a pH of 5 to a pH of 3, the hydronium ion concentration is			2) red with litmus	nus			
	(1) 0.01 of the original content			3) colorless with ph	nenolphthalein		
	(2) U. For the original content		0	pink with phenol	phthalein - Tul	LNS PINK	
	(3) 10 times the original (4) 100 times the original				ı~ A	BASE	
	2 41 5				•		
	2 7 3						
	10× /0X	= 100 X					

25. The results of testing a colorless solution with three indicators are shown in the table below.

A BASE

Indicator	Result		
red litmus	blue		
blue litmus	blue		
phenolphthalein	pink		

## A METAL + OH = BASE

Which formula could represent the solution tested?

- 1) NaOH(aq)
- (3)  $C_6H_{12}O_6(aq)$
- 2) HCl(aq)
- (4)  $C_{12}H_{22}O_{11}(aq)$
- 26. Which equation represents a neutralization reaction?
  - ((1))NaOH + HCl → NaCl + H<sub>2</sub>O
  - (2) 2Na + 2H<sub>2</sub>O → 2NaOH + H<sub>2</sub>

TABLE I

- (3)  $Zn + CuSO_4 \rightarrow ZnSO_4 + Cu$
- (4) AgNO<sub>3</sub> + NaCl → AgCl + NaNO<sub>3</sub>
- 27. Which products are formed when an acid reacts with a base?
  - (1) an alcohol and carbon dioxide
  - (2) an ester and water
  - (3) a soap and glycerine
  - (4) a salt and water

TABLET

- 28. Which equation represents a neutralization reaction?
  - (1) CaO + H<sub>2</sub>O → Ca(OH)<sub>2</sub>
  - (2)  $2HCI + Zn \rightarrow ZnCI_2 + H_2$
  - (3)  $H_2SO_4 + CaCO_3 \rightarrow CaSO_4 + H_2O + CO_2$
  - (4)  $HNO_3 + KOH \rightarrow KNO_3 + H_2O$

AUD BASE SALT

- 29. What is the pH of a solution that results from the complete neutralization of an HCl solution with a KOH solution?
  - (1) 1

- (3) 10
- 14,0
- (4) 4
- 30. How many moles of KOH are needed to exactly neutralize 500. ml of a 1.0 M HCl solution?
  - (1) 1.0

(3) 0.25

(2) 2.0

- 31. How many milliliters of 0.2 M NaOH are required to exactly neutralize 40 milliliters of 0.1 M HCI?
  - (1) 10

(3) 40

20


- (4) 80
- 32. A 9.0 milliliter sample of HCl solution is exactly neutralized by 6.0 milliliters of a 3.0 M NaOH solution. What is the concentration of the HCl solution?
  - (1) 1.0 M

(3) 1.8 M

(2) 2.0 M


(4) 1.5 M

33. The diagram below represents the menisdus on an acid and a base buret at the endpoint of a titration in which 0.10 M NaOH was used to neutralize an unknown concentration of HCI.



If the solution level in each buret was 0.00 milliliter at the start of the titration, what is the molarity of the unknown HCI solution?

- (1) 1.2 M
- (3) 0.30 M
- (2) 0.13 M
- (4) 0.090 M
- 34. The diagram below shows NaOH(aq) being added to HCI(aq). A few drops of phenolphthalein were added to the flask before the titration was started.



The endpoint in this titration is reached when the solution in the flask appears

(1)) pink

(2) colorless

(3) blue (4) yellow



- 35. If 20. milliliters of 4.0 M NaOH is exactly neutralized by 20. milliliters of HCI, the molarity of the HCI is
  - (1) 1.0 M

(3) 5.0 M

(2) 2.0 M

(4)) 4.0 M




TABLE M

QUESTIONS 30 - 35

Ma Va = My Vb

MOLARITY X VOLUME = MOLARITY X VOLUME
ACED BASE

30. HC1 = KOH

ma va mb Vb

/ m × 500 ml = m × 500 ml

500 m = 1 moverity

M = ? moles molary

. 5 MOLES OF KOH

31. HC1 = NaOH

Ma Va = Mb Vb

. Im · 40me = 2m · V

4 = .2 V

4/.2 = V

20 = V

32, HC1 = NaOH

m.9m1 = 3m.6m1

9m = 18

m = 20

33. It c1 = NaOH

ma Va Mb Vb

m · .9m1 = .1m · 1.

.9m= .12

m = .13m

34. TABLE M

35. HU = NaOH

ma Va mb Vb

M. 20 ml = 4/mx 20 ml

20m = 80

(m = 4 m)